






Get the most out of your qPCR data with the market leading GenEx software

GenEx, the popular software for the processing and analysis of real-time qPCR data, provides a multitude of functionalities for the qPCR community. Its capabilities range from basic data pre-processing and management to advanced cutting-edge data analysis. GenEx provides methods to select and validate reference genes, classify samples, group genes, monitor time dependent processes, and much more.




Comprehensive workflows for NGS: From QC to figures in one go

- Web-based: No software to install
- Easy to use: Run analyses in minutes
- Reports: Understand your data with a quick glance

You register and purchase analysis time on: tataa.basepairtech.com

Contact us for a free trial: services@tataa.com



HOME FEATURES PLANS AND PRICING BLOG CONTACT US

Analyze your NGS data in hours

Upload your raw data, choose among several high quality workflows for RNA-Seq, ChIP-Seq, DNA-Seq, etc. and get results in an hour.

[Start free trial](#) [Learn more](#)

Trusted by scientists at



News

[Biomine "datautvinning för identifiering, urval och validering av biomarkörer" beviljas 12 miljoner av KK-stiftelsen. Biomine koordineras av Högskolan i Skövde och TATAA Biocenter är en av deltagarna.](#)

[CAMP – Centre for Advanced Medical Products, beviljas 48 miljoner av Vinnova. CAMP koordineras av Umeå universitet och TATAA Biocenter är en av deltagarna.](#)

Don't forget our summer offers!

[19% discount* on our new TATAA GenMaster Mix – your best choice for routine qPCR](#)
[> TATAA SYBR GenMaster Mix – The easy-to-use qPCR fastmix for all your gene expression experiments](#)

[> TATAA Probe GenMaster Mix – The easy-to-use qPCR fastmix for all your gene expression and genotyping experiments](#)

[50% summer discount** on TATAA GrandPerformance Cancer Assays](#)

TATAA GrandPerformance assays are the only on the market validated for single cell profiling. Assays from the TATAA Cancer panel can be mixed in any combination and used

in a preamplification workflow offering a sensitivity down to a few transcripts under optimum conditions.

[Buy GrandMaster and receive assay design and validation for 150 €!](#)

At TATAA we celebrate Midsummer by offering all customers that order a [TATAA GrandMaster Mix](#) assay design and validation from only 150 €* through our partner!

Let us help you with your nucleic acid analysis!

Find high quality products for your qPCR work in our [webshop](#)

To get expert advice and consultation, you are welcome to use our [commissioned services](#)

Get more knowledge and help yourself through our [hands-on courses](#)



TATAA BIOCENTER AB
Odinsgatan 28
SE-411 03 Gothenburg